

RCC-MS-GMD23G-01R

23.8" FHD Medical Monitor For General Clinical Usage



Product Features

- IP65 in front & IPx1 whole unit
- 20 meter DVI & HDMI long cable solution
- Calibration mode (Optional)
- DICOM compliant (Optional)

Specification

Audio

Speaker Output 2 x 0.5W Speaker

Certification

Certification Information FCC Class B, VCCI, CCC
IEC/EN/BS EN60601-1, IEC/EN/BS EN60601-2
cTUVus certified (ANSI/AAMI ES60601-1, CSA C22.2 NO.60601-1)
DIN6868-157 compliant (optional)

Panel

LCD Size 23.8"

Display Type FHD

Resolution 1920 x 1080

Luminance AG glass 230 nits
AG PCAP 220 nits (optional)

Contrast ratio 3000: 1

Viewing angle 178°/ 178°

Response time 16 ms

Backlight LED

Touch Type	Projective capacitive multi-touch up to 10 points
------------	---

Display color	16.7M
---------------	-------

DICOM compliant	Barten curve (optional)
-----------------	-------------------------

I/O

Input Interface	1 x VGA 1 x DVI 1 x HDMI 1 x Display port
-----------------	--

User control and indicator	1 x Power indicator & Sensor keys multi-language OSD: contrast, brightness, auto adjust, clock, phase, color temp, hor position, ver position, input select , key lock
----------------------------	--

Mechanical & Environmental

Weight	6.8 kg
--------	--------

IP Rating	IP65 in front & IPx1 whole unit
-----------	---------------------------------

Mounting	VESA 100 x 100
----------	----------------

Construction	Plastic
--------------	---------

Power Requirement	AC/DC Adapter Input: 100~250Vac, 0.6A, 60/50Hz universal Output: 24V, 2.71A
-------------------	---

Power Input	24V DC
-------------	--------

Dimension (mm)	587 x 367.5 x 59.5
----------------	--------------------

Power consumption	On:<40W; active off: <5W
-------------------	--------------------------

Power management	EPA Energy Star compliant
------------------	---------------------------

Color Mode

Color temperature	Cool/neutral/warm/user
-------------------	------------------------

gamma	1.8/2.0/2.2/2.4/linear
-------	------------------------

Calibration mode	Communication protocol is ready to calibrate the monitor through DDC CI (optional)
------------------	--

Ordering Information & Packing List

RCC-MS-GMD23G- 01R	Medical Monitor with AG Glass, 23.8"/LED/AG Glass/Plastic Case/1920x1080/230nits/3000 Contrast/ 0.5W Speakerx2/DVI/DP/HDMI/VGA Packing List: Medical Monitorx1, Power Adapterx1, HDMI 1.8M Cables1
--------------------	--

RCC-MS-GMD23T- 01R	Medical Monitor with AG PACP, 23.8"/LED/AG PCAP/Plastic Case/1920x1080/220nits/3000 Contrast/ 0.5W Speakerx2/DVI/DP/HDMI/VGA Packing List: Medical Monitorx1, Power Adapterx1, HDMI 1.8M Cables1, USB Touch Cables1
--------------------	---

Optional Items

Optional Items	AG PCAP touch
	DICOM
	Calibration mode
