

HID-2100

21.5" FHD Medical Monitor For General Clinical Usage



Product Features

- Front IP65 protection
- 250 nits brightness
- HDMI & DP video input

Specification

Audio

Speaker Output 2 x 2W Speaker

Certification

Certification Information CE, FCC

Panel

LCD Size 21.5"

Display Type FHD

Resolution 1920 x 1080

Luminance 250 cd/m²

Contrast ratio 1000: 1

Viewing angle 178°/ 178°

Response time 22 ms

Backlight LED

Touch Type	Projective capacitive multi-touch up to 10 points
-------------------	---

I/O

Input Interface	1 x VGA 1 x HDMI 1 x Display port
------------------------	---

Mechanical & Environmental

Weight	5.5 kg
IP Rating	IP65 in front
Mounting	VESA 75 x 75, 100 x 100
Construction	Plastic
Power Requirement	AC/DC 120W Medical Adapter
Power Input	24V DC
Dimension (mm)	539.6 x 342.6 x 45.2

Ordering Information & Packing List

HID-2100-252-A1R	21.5 PCAP Multi-touch Monitor 1920x1080/DP/HDMI/VGA/Line in/Line out/DC out 24V 30W/2W speaker Front OSD key/2xUSB/USB type B for Touch/ 24V DC in/120W Medical Power adapter
-------------------------	--
